

Technical product specification

Product name	semperguard Style	Version / Index no:
Spec code	NOF-035BK-N-3CZ	semperguard Style_NOF-035BK-N- 3CZ Version A November 2022 EN
Date of issue	24.11.2022	

General information

Type single use examination and disposable protective glove, non sterile

Labelling information printed on dispenser box
Shape ambidextrous - straight fingers

Material Nitrile Butadiene Rubber (NBR) [not made with natural rubber latex]

ColourblackInsidepowder freeOutsideno treatment

Cuff / surface rolled cuff / finger textured

Shelf life 3 years

Available sizes XS (5-6) S (6-7) M (7-8) L (8-9) XL (9-10)

Dimensions, physical properties and biocompatibility

Glove length median ≥ 240 mm (according to EN 455-2)

Minimum wall
thicknessat finger0.14 mm (double measured) / 0.07 mm (single measured)at palm
at cuff0.12 mm (double measured) / 0.06 mm (single measured)0.10 mm (double measured) / 0.05 mm (single measured)

Glove width according to EN 455-2: S 80 ± 10 mm, M 95 ± 10 mm, L 110 ± 10 mm, XL ≥ 110 mm

Force at Break median ≥ 6 N (during shelf life according to EN 455-2)

Tensile Strength min. 14 MPa after aging (according to ASTM D6319)

Elongation at Break min. 400% after aging (according to ASTM D6319)

Residual powder / ≤ 2 mg (according to EN 455-3)

Powder content

Performance requirements and inspection levels

Freedom from holes (Barrier) AQL ≤ 1.5

(as per EN 455-1, sampling in accordance with ISO 2859-1, G-1)

Dimensions and physical properties AQL 4.0

(as per ASTM D6319, sampling in accordance with ISO 2859-1, S-2)

Standards, guidelines & quality certificates

Quality certification ISO 9001, ISO 13485, ISO 14001, ISO 45001

Conformity to regulations - Medical Device Regulation (EU) 2017/745: Class I

- PPE Regulation (EU) 2016/425: Category III

- Regulation (EC) 1935/2004 on Food Contact Materials

Conformity to standards EN 420 / EN ISO 21420, EN ISO 374-1 (Type B), EN 374-2, EN 16523-1, EN

374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1,

ASTM D6319, ASTM F1671/F1671M



Technical product specification

Product name	semperguard Style	Version / Index no:
Spec code	NOF-035BK-N-3CZ	semperguard Style_NOF-035BK-N- 3CZ Version A November 2022 EN
Date of issue	24.11.2022	

Instructions and additional statements

Storage instruction Store in original packaging in a dry and dark place at 10 °C to 30 °C. Refer to

guidelines of storage of rubber products as described in ISO 2230:2002. Ensure that storage area is kept cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper ions discolour the glove. Protect gloves against ultraviolet light sources, such as sunlight and oxidizing agents. Storage above 30 °C will lead to accelerated aging and should be

avoided.

Cautionary statement and ingredient information

This product contains accelerators (Dithiocarbamate types, Zincmercaptobenzothiazol) not to be used in a hypersensitivity of these

substances.

For further information, a list of substances contained in the glove is available

upon request.

Reporting system

Medical device vigilance and reporting system

According to the official reporting criteria of the Medical Device Regulation, incidents caused by examination gloves must be reported immediately to our

M. Schirmbrandt

Product Management

Medical Device Reporting team. E-Mail:

Head of Quality Management and

sempermed.complaints@semperitgroup.com or Tel.: +43 2630 310 0

A. Wöss Director

Remark

Replaces all previous versions.
All standards references refer to the date of document issue.

Regulatory Affairs

B. Sebauer