

sempermed® 
derma PF

DERMA PF FOR SENSITIVE HANDS

HIGHLY SKIN-FRIENDLY – the Semperit Post Production Process
EASY TO DON – patented synthetic inner coating
MAXIMUM COMFORT – fatigue-proof working



FOR SENSITIVE HANDS

The Sempermed® Derma PF offers good reasons to switch to powder-free gloves.

■ HIGHLY SKIN-FRIENDLY

This powder-free glove provides improved safety for staff, patients and visitors as there are no powder-related dermatological, allergic and toxic health risks.

■ EASY TO DON

With its synthetic inner coating the glove is easy to put on due to its patented net structure.

■ SURE GRIP

The Sempermed® Derma PF has more grip and thus combines safety, tactical properties and instrument guidance.

■ GOOD VALUE FOR MONEY

The long-term cost reduction is significant, if you take the very high follow-up costs of powder damage into consideration. The cost of purchasing Sempermed® Derma PF are minimal by comparison.



TECHNICAL DATA

sempermed® derma PF

TYPE	Sterile surgical glove for single use, powder-free, with synthetic lining		
MATERIAL	Natural rubber latex (NR)		
COLOR	White		
GLOVE SHAPE	Fully anatomical with rolled rim		
SIZE/OVERALL LENGTH as per EN 455-2	5 ½, 6 and 6 ½ 7, 7 ½ and 8 8 ½ and 9	min. 265 mm min. 275 mm min. 280 mm	
WALL THICKNESS measured in single layer in the palm area	0.16 - 0.20 mm		
IMPERMEABILITY / FORCE AT BREAK as per EN 455-1 / as per EN 455-2	AQL 1.0 / ≥12 N		
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years		
STERILISATION	Radiation STERILE R		
LABELLING	EN 1041, EN 420, ISO 15223		
PURPOSE as set out in MDD 93/42/EEC as set out in PPE 89/686 EEC	Medical device class IIa Single-use protective glove – PPE category III (protection against chemical substances for limited time) ¹⁾		
PACKAGING	Left and right glove with turned up cuff in fibre-free inner pouch, ozone-tight, sealed in flat medical peel pack. In dispenser carton with sterilisation indicator: 50 pairs In transport carton with sterilisation indicator: 300 pairs		
ARTICLE NUMBERS	Size 5 ½ 822751522 Size 6 822751602 Size 6 ½ 822751622 Size 7 822751702	Size 7 ½ 822751722 Size 8 822751802 Size 8 ½ 822751822 Size 9 822751902	



CONTACT US!

Semperit Investments Asia Pte Ltd
8 Jurong Town Hall Road,
#29-03 to 06 JTC Summit, Singapore 609434
e-mail: sempermed@semperitgroup.com
www.sempermed.com

Semperit Technische Produkte Gesellschaft m. b. H.
Segment Sempermed
Austria · 1030 Vienna · Modecenterstraße 22
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630
e-mail: sempermed@semperitgroup.com · www.sempermed.com

¹⁾ Please note the Sempermed Chemical Resistance list, available on www.sempermed.com/sempermed_informs/chemical_resistance

IMPORTANT NOTE: Actual product information is available under www.sempermed.com. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. Semperit does not accept any liability for incorrect use of the gloves. In case of doubt, obtain expert advice before use. The information and classification correspond to the latest status prior to printing. Technical details are average values from production and may vary in individual cases. Subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock.**