

### PRODUCT SPECIFICATION

PROTECTS HYGIENIC LATEX powdered

#### GENERAL INFORMATION

**Type:** single-use disposable protective glove  
**Marking:** all information on dispenser box  
**Shape:** ambidextrous – straight fingers  
**Material:** natural rubber latex (NRL), powdered  
**Colour:** natural white  
**Cuff/ Surface:** rolled cuff/ finger textured  
**Thickness:** (double measured)



Finger: min. 0.16 mm (AQL 4.0)  
 Palm: min. 0.16 mm (AQL 4.0)

**Length:** S median 220 mm  
 M, L, XL median 230 mm

**Storage:** Protect the gloves from heat, humidity, strong light and ozone. Storing at 10°C – 30°C in a dark and dry place with 60% - 80% humidity is recommended. Date of expiry on dispenser box and transport carton.

**Article numbers:** S (6-7) 8139 95823  
 M (7-8) 8139 95825  
 L (8-9) 8139 95827  
 XL (9-10) 8139 95829

#### STANDARDS, GUIDELINES AND QUALITY CERTIFICATION

**Quality certification:** ISO 9001  
**Conformity to guidelines:** Directive 89/686/EEC on personal protective equipment  
 CE Category I for minimal risks

**Conformity to standards:** EN 420 EN 374-2: specified  
 level 1: AQL 4.0 (S4)

**Conformity to food regulation:** 1935/2004 

#### WARNING AND REPORTING SYSTEM



**Potentially allergenic ingredients:** Dithiocarbamate type, Zinc-mercaptobenzothiazol  
 (For further information the actual version of the ingredients list is given on request.)  
 Caution! This product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses!

#### PACKAGING


**Packaging volume:** Dispenser box S, M, L 200 pcs  
 Dispenser box XL 180 pcs  
 Dispenser boxes per transport carton 10

**Packaging dimensions:**

| in mm            | width | length | height |
|------------------|-------|--------|--------|
| Dispenser box    | 128   | 245    | 80     |
| Transport carton | 255   | 418    | 264    |

| Established (PM/DM2):   | Reviewed (QA/DM):   | Approved (HD/DM):  | Customer approval |
|---|---|--|-------------------|
|  |  |  |                   |
| Mag. M. Schirmbrandt  | DI (FH) B. Sebauer  | DI C. Eichler  |                   |

PROTECTS\_HYGIENIC\_LATEX\_PP\_EN\_Version A, 2012-05-25, replaces all previous versions

Semperit Technische Produkte Gesellschaft m.b.H. • Modecenterstrasse 22 • A - 1031 Wien,  
 t+ 43 / 1 / 79 777 – 520, Fax: + 43 / 1 / 79 777 – 630  sempermed@semperit.at; www.sempermed.com