



150250

Tork Couch Roll Advanced

Protect your patient couch with the disposable Tork Couch Roll Advanced. The perforated couch roll is ideal to maintain a high level of hygiene in a healthcare environment. Rolls can be used in Tork Couch Roll Dispenser.

- Perforation for easy tear off and consumption control
- Strong and easy to roll, for optimal protection
- Better resistance and higher comfort due to 2-ply

Product certificates



The mark of
responsible forestry



Think ahead.

Product and delivery data

System C1	Quality Advanced	Roll length 54.5 m
Roll width 50 cm	Roll diameter 12 cm	Number of sheets 165
Sheet length 33 cm	Core inner diameter 5.1 cm	Ply 2
Number of sheets 165	Embossing Yes	Colour White

Delivery data

	Consumer units (CON)	Transport units (TRP)	Pallet (PAL)
EAN	7310791230351	7310791230344	7322540904192
Packaging material	Plastic	Plastic	
Pieces	1	9 (9 CON)	162 (18 TRP)
Transport units per layer			9
Layer per pallet			2
Gross weight	1023.1 g	9.3 kg	167.09 kg
Height	500 mm	500 mm	1161 mm
Length	120 mm	360 mm	1200 mm
Width	120 mm	360 mm	1000 mm
Net weight	980.1 g	8.8 kg	158.78 kg
Volume	7.2 dm ³	64.8 dm ³	1.17 m ³

Compatible products



[129184](#)
Tork Couch Roll
Dispenser

[129183](#)
Tork Couch Roll
Dispenser